

Field Change Order

FCO 70400042

Optimus RAD release upgrade to 3.8 USA only



BuckyDiagnost & Cosmos

DOCUMENT HISTORY:

Revision	Revision date	Reason of changes	
AA	Jan. 2011	First issue	4512 980 65771 REV AA
AB	April 2011	Revised document	4512 980 65772 REV AA

CSIP Security Labeling: CSIP Level 1:

This document and the information contained in it is strictly reserved for current Philips Healthcare ("Philips") personnel, Philips licensed representatives and Philips customers who have purchased a valid service agreement for use by the customer's designated in-house service employee on equipment located at the customer's designated site. Use of this document by unauthorized persons is strictly prohibited. This document must be returned to Philips when the user is no longer licensed and in any event upon Philips' first written request.

Proprietary Notice:

This document and the information contained in it is proprietary and confidential information of Philips Healthcare ("Philips") and may not be reproduced, copied in whole or in part, adapted, modified, disclosed to others, or disseminated without the prior written permission of the Philips Legal Department. Use of this document and the information contained in it is strictly reserved for current Philips personnel and Philips customers who have a current and valid license from Philips for use by the customer's designated in-house service employee on equipment located at the customer's designated site. Use of this document by unauthorized persons is strictly prohibited. Report violation of these requirements to the Philips Legal Department. This document must be returned to Philips when the user is no longer licensed and in any event upon Philips' first written request.

Warranty Disclaimer Language:

Philips provides this DOCUMENT without warranty of any kind, implied or expressed, including, but not limited to, the implied warranties of merchantability and fitness for a particular purpose.

Limitations of Liability Language:

Philips has taken care to ensure the accuracy of this document. However, Philips assumes no liability for errors or omissions and reserves the right to make changes without further notice to any products herein to improve reliability, function, or design. Philips may make improvements or changes in the product(s) or program(s) described in this document at any time.

APPLIES TO:

Geography:	USA only
Traceable item identification (12NC):	9890 000 02001
Range of serial numbers:	SN Optimus generator: ≥ 970218 and ≤ SN061263 (delivered: from May 1997 - January 2007)
Physical sub, main block(s) and system code(s) identification where the affected item could be present:	70410/11/12/14/15, 704030/35, 704020/21/22, 704031/32, 704060/62, 704016/17/18 Main block: PB000135 Sub block: PB010005, PB010032, PB010033
Any other type of system identification number used by the BG/BU:	Optimus generator with FW version ≤ 3.6 (delivered: from May 1997 - January 2007) Control panel 9890 000 02403/-04/-05/-06/-07/-08 (delivered with SN: ≥ 9701016 and ≤ SN070056)

REMOTE SOFTWARE INSTALLATION:

Is this a software "ONLY" FCO: Yes

Can the software be installed on the equipment remotely: No

Procedure: Update of the SW with latest release

OPPORTUNITIES FOR FCO IMPLEMENTATION EFFICIENCY:

N/A

SUPERSEDES INFORMATION:

N/A

INTRODUCTION:

Symptom : Exposure parameters like kV, mA, mAs or ms seem to change without a button having been touched.

List FPR(s) solved (if applicable)
PR # 333162; 342432; 366270

Cause : If one pushes continuously (longer than a second) the + or – button to change any of the exposure parameters, the generator SW switches into an auto stepping mode, so that the data continues to change (up or down). This function is specified to stop whenever you release the button.

If you are in the auto stepping mode (or within less than 200ms after release of this function) AND also the preparation switch is pressed, the auto stepping mode does NOT stop.
This might lead to unwanted parameter settings.

Remedy : Upgrade of the generator SW to latest version release 3.8

MANPOWER / TIME TO COMPLETE:

1 service engineer / 1.0 hour

APPLICATION TRAINING REQUIREMENTS:

No

COMPLIANCE TESTING:

No

TOOLS & TEST EQUIPMENT:

Tools	Tool code / 12NC
Standard tool kit	TC129
Service engineer PC, IST and AGenT latest version	TC 092

MODIFICATION KIT / PARTS REQUIRED:

Modification kit and/or parts 12NCs	Description
N/A	N/A

New firmware included in: AGenT 5.3.1, use Zeppelin tool CD or download from InCenter.

[Link to Zeppelin ToolBox](#)

FCO KIT CONTENTS:

Non-traceable items

Item 12NC	Description
N/A	N/A

Non-serialized trace items

System code	FCO document	Outbound item	Inbound item	Parent material number	Software release description (N/A when HW)	Software patch level (N/A when HW)
Trace items (OIs) that must be installed/exchanged/removed (MP1 OI-table)						
N/A	N/A	N/A	N/A	N/A	N/A	N/A

Ordering information: Order the indicated material/kit according to the standard local service logistic procedures.

INSTALLED BASE REGISTRATION:

This FCO has impact on the installed base registration.

Make sure your local installed base registration is or gets updated with the 12NC numbers (and/or serial numbers) as stated in the **trace item tables** under "FCO Kit contents".

PROCEDURE:

1 PREPARATIONS

1.1 Off-site preparations

- Check the system status.
 - System compatibility
 - Other applicable FCOs
- Arrange the visit.
- Arrange the required tools and required firmware.

1.2 Back-up

NOTICE



Use proper ESD grounding techniques when handling components
Wear an antistatic wrist strap and use an ESD-protected mat.
Store ESD-sensitive components in antistatic bags before placing them on any surface.

1.3 Preparing the generator

1.3.1 Preparing generators without a CAN interface

The loading process can be started once relay ENK1 has been energized. Proceed to 1.3.3.

1.3.2 Preparing generators which are connected via a CAN interface

BuckyDiagnost TH and TH2

- Switch OFF the generator.
- Disconnect the following plugs:

System	Connector		
	EZX23 signal bus	EZX42 or EZX42-1 system CAN	EZX43 or EZX43-1 system CAN
BuckyDiagnost TH / TH2	X		X

Important



The download procedure must not be started before relay ENK1 has been energized at least two minutes after the generator has been switched ON.

1.3.3 Establishing the PC-generator connection

- If IST is not already started on the PC, start it now.
- Switch the generator ON.
The download process can be started once relay ENK1 has been energized.
- Establish connection between PC and generator with the data cable.
- Start AGenT by clicking on the respective icon.



Tip

CU backup is the battery buffered configuration excluding tube adaptation data. Firmware release 3.x is loaded and memorized in flash PROMs. Loading the flash firmware does not influence the CU backup configuration.


1.3.4 Saving CMOS (CU backup), APR and tube load statistics



Tip


The backup is a safety backup which shall replace older ones which might already be obsolete from the configuration and in case it might be needed for a CU fault. A new backup after the upgrade procedure is not necessary and it is not necessary to restore it after the firmware has been updated.

- Check the actual level of release 3.
- Select menu:
AGenT / Fault Find / Power ON Results / SW/HW Versions.




R3V1L8 = Release 3 Version 1 **Level 8**


- Select menu:
Agent / Acceptance / Backup.



- Download the actual CMOS (filename e.g. CU97xxxx. TDL = serial number of generator)
It takes about 15...20 minutes to save the data to the disk.
The default backup name: **CUBACKUP.TDL** can be changed into any other file name.
The path (hard disk) is automatically taken into account.
- At the end of the downloading process wait until the following message appears.



- Reset the generator.



- Click Yes to restart the connection if the actual release is 3.5 or 3.6 and proceed with **2. FCO Implementation**.
- Click No if the actual release is < 3.5.

The maximum value for density correction has been limited from Optimus RAD/RF release 3.5 onwards.

At Bucky auxiliaries the sum of all density corrections is

max. +/- 2 steps of 25% (+/- 4 formerly)
or
+/- 4 steps of 12% (+/- 8 formerly)
or
+/- 8 steps of 6% (+/- 16 formerly).

This is valid for all auxiliary settings, individual APR settings and patient size corrections.


This requires that the APR is downloaded on generators with release version <3.5 to have the data available for maintenance in case there are higher correction values programmed than possible after the upgrade.

Higher correction values are automatically limited after the flash load process.

- Use program APRMAN rel. 2.1 from the Zeppelin platform to download APR.
- Carry out the required modifications for the + density corrections > 2x25%. Increase the background organ mAs and select a slower film speed. The file can later be reloaded to the generator after the flash upgrade to release 3.8.
- Stop the APRMAN program.

From release 3.5 onwards the tube load counter has more detailed information. The content of the counter of the predecessor releases gets lost with the 3.8 update.

- Start AGenT.
- Select menu:
AGenT / Acceptance / Inspect / Tube 1 (2, 3).
- Record the entries in the system logbook.




- Keep AGenT connected for the next step.


2 FCO IMPLEMENTATION

2.1 LOADING THE GENERATOR CU FIRMWARE RELEASE 3.8

- Select menu:
AGenT / Program / Update Generator Firmware (XRG 90 RAD/RF,C)




- Select the respective update file **OMCUR3L8.TDL** and click on "Open" with the left mouse button.



The reset can be performed within the next forty seconds, either with PCB EZ 139 S1 or with the ON button of the control module.

During the update process a progress bar is displayed on the screen which indicates how much of the update is completed.



Depending on the type of PC; data transmission takes 15 ... 30 minutes, in some cases about an hour. During this process all red LEDs of the function units are blinking, the CU LED is permanently on.

Tip




When the data transmission to the generator is completed, the message to wait for two minutes appears on the screen.



Important




This process must under no circumstances be disturbed!
At the end of this sensible procedure "Flash loaded ok" appears on the screen.
Only now the AGenT program can be terminated Text.



Typically the function units stop blinking, CU blinks. Once CU stops blinking as well the generator should come up normally.

- Select menu:
Agent / Faultfind / Power ON Results / SW/HW Versions



to confirm that release 3.8 has successfully been loaded.

2.2 FINAL WORK AND REMARKS ON THE CHANGES

Generators which are connected via a CAN interface:

- Switch OFF the generator.
- Reestablish the signal bus connector EZX23 and the CAN connector EZX43/43-1.

There are no changes which can be recognized at first sight.

In the error log index (path: *AGenT / Fault Find / Error log*) there is a 00XQ entry if the original firmware was <3.5 before.

It indicates "Tube statistic data invalid" as the tube load table changed its format.

Release 3.8 has a detailed table now.

All tube counters of the predecessor release are empty after flash load of release 3.8 if the original firmware was <3.5 before.

Details and a table of the new counter table can be found in chapter 3. "Exposure counter" at the end of this manual.

2.2.1 Settings

No programming screen has to be adapted to the new level 8 of release 3 after the flash load, but some changes of settings have to be carried out:

Tomography auxiliary settings have to be modified for CAN controlled systems such as

Bucky TH, TH2:

RGDV 2 Data Set A

(see pages 2Z-2.2 + 2Z-2.5 + 2Z-2.8 + 2Z-2.9)
at the end of this document

Exposure series / Tomo movement: **No**

RGDV 2 Data Set B

(see pages 2Z-2.2 + 2Z-2.5 + 2Z-2.8 + 2Z-2.9)
at the end of this document

Underexposure display (non automatic techniques): **Yes**

2.2.1.1 Names in the screen

AGenT / Program / RGDV set A + B / RGDV 1 ... 8 / Data Set A
of the programmed "Mounted radiographic controller" changed from
Bucky Controller 1 into Bucky Ctrl. 1/Dig.Diag.

2.2.2 APR data

- If APR data have to be adapted to the different maximum dose corrections, re-load the modified APR data file to the generator now.

2.2.3 Error log

The error log index can be erased if there are no entries which require service work:

AGenT / Fault Find / Error Log

- Click on "Clear" with the left mouse button.

3 EXPOSURE COUNTER

- Before handing over the generator to the customer, read the exposure counter.

Use menu:

Acceptance / Inspect / Tube 1 ... 3 / Type and statistic of Tube 1 ... 3

Record the figure in the table below.



Tip

Tube load statistics variables written on a grey background and marked by a "*" are visible but do not affect the functions of this generator RAD type.
(They are made for generators R/F version).

Tube load statistic variable	Unit	Tube1	Tube 2	Tube 3
Reset date	dd.mm.yy			
Last update	dd.mm.yy			
Preparation time large focus	s			
Preparation time small focus	S			
Preparation time vario focus	S			
Preparation counter large focus	1			
Preparation counter small focus	1			
Preparation counter vario focus	1			
* Fluoro time	min			
* Fluoro counter	1			
Exposure counter large focus	1			
Exposure counter small focus	1			
Exposure counter vario focus	1			
Overload exposures counter large focus	1			
Overload exposures counter small focus	1			
Overload exposures counter vario focus	1			

- The tables should be reset whenever the tubes are being replaced.
Use menu:
Acceptance / Inspect / Tube 1 ... 3 / Type and statistic of Tube 1 ... 3
- Click on "Reset" with the left mouse button.
- Record the figure in the table above.

Type and Statistic of Tube 1 [?] [X]

Menu/Acceptance/Inspect/Tube 1

[OK] [Save] [Reset] [Help On]

Type of Tube

SR033100-ROT350

Tube Statistic

reset date	30.05.01
last update	26.03.04
preparation time large focus	16390
preparation time small focus	6525
preparation time vario focus	2
preparation counter large focus	4094
preparation counter small focus	5205
preparation counter vario focus	1
fluoro time	85
fluoro counter	548
exposure counter large focus	5047
exposure counter small focus	5364
exposure counter vario focus	1
overload exposure counter large focus	0
overload exposure counter small focus	4
overload exposure counter vario focus	0

Explanation:

Reset date / Last update:

Reset of date and date of last update of the tube statistics.

Preparation time:

The sum of all preparation times per focus.

Preparation counter:

Count of the occurrences of transition STANDBY or FLUORO to PREPARATION per focus.

* Fluoro time:

The sum of all fluoro times.

* Fluoro counter:

Count of the fluoro commands.

Exposure counter:

Count of the exposures per focus (including the overload exposures).

Overload exposures counter:

Count of the exposures at overload conditions of the tube.

4 FINISHING WORK

- Close the cabinet.

INSTALLED BASE REGISTRATION – FSE NOTE:

This FCO has impact on the installed base registration so make sure you update your local installed base registration.

PARTS DISPOSAL:

N/A

DOCUMENTATION:

- Log this action in the section "History Record" of the System Reference Manual.
- File this FCO in the section "Service Information" of the System Reference Manual.
- Fill out the attached Action Notification Report and send it to your SSD Customer Support Manager.

ACTION NOTIFICATION REPORT:

- If required, fill out the attached ANR (Action Notification Report) and send it to your local GS&S Key Market / Country Customer Services or FCO manager.

Bucky Diagnost TH with Bucky-Controller		Name :		Bucky		Tomo		Bucky wall stand		Free cassette		D 76 Scopo		D 76 BV - DSI			
- D76 : Exposure Scopo / BV-DSI		Desk :		RGDV1		RGDV2		RGDV3		RGDV4		RGDV5		RGDV6		RGDV7	
Data Set A :				1		1		1		1		1		1			
- Room :																	
- Tube :				1		1		1		1		2		2			
- Release circuit number :				do not care		do not care		do not care		do not care		1		2			
- Enable handswitch at generator desk :				yes		yes		yes		yes		no		no			
- Synchronizer present :				yes		yes		yes		yes		yes		yes			
- Exposure switch type :				double step		double step		double step		double step		double step		double step			
- Bucky format density correction :				0		0		0		0		0		0			
- Cone density correction :				0		0		0		0		0		0			
- Dose measurement input :				EZ X21		none / [EZ X21]		EZ X31		none		EZ X22		EZ X41			
- Dose measurement sensor :				Bucky amplimat		Bucky amplimat		Bucky amplimat		(Bucky amplimat)		Scopo amplimat		photo sensor / amplimat input			
- Exposure series/Tomo movement :				no		no		no		no		yes		yes			
- Release delay :				enable		enable		enable		enable		enable		enable			
- Mounted radiographical controller :				Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		none		none			
- Release circuit adaptation unit :				none		none		none		none		1WB		1WB			
- Mounted tomo extension :				0		0		0		0		0		0			
- Medium II format kV correction (dose equiv. steps) :				0		0		0		0		0		0			
- Medium II format density correction (6% steps) :				0		0		0		0		0		0			
- Medium II format mAs correction (6% steps) :				0		0		0		0		0		0			
- Small II format kV correction (dose equiv. steps) :				0		0		0		0		0		0			
- Small II format density correction (6% steps) :				0		0		0		0		0		0			
- Small II format mAs correction (6% steps) :				0		0		0		0		0		0			
Data Set B :																	
- Used for tomo :				no		yes		no		no		no		no			
- Used for fluoroscopy :				no		no		no		no		yes		yes			
- CT add on :				no		no		no		no		no		no			
- Disable time override :				no		no		no		no		no		no			
- Tube power factor :				100 %		100 %		100 %		100 %		100 %		100 %			
- kV steps :				Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)			
- mAs steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- mA steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- time steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- Density steps :				12 % 1)		12 % 1)		12 % 1)		12 %		12 % 1)		12 %			
- Density correction (6% steps) :				0		0		0		0		0		0			
- Underexposure display :				yes		yes		yes		yes		yes		yes			
- Tube overload protection :				on		on		on		on		on		on			
Bucky / Scopo 1WB / Decade Bucky 1 (WBX11) :				Bucky / Tomo 1WA : Decade Bucky 1/2		WAX41		WAX41		WAX42		Bucky / Tomo 1WA : Tomo time [s] :					
				Tomo mode switch								Tomo time 1 :		Tomo time 5 :			
				Bucky RGDV - switch related								Tomo time 2 :		Tomo time 6 :			
				Bucky RGDV								Tomo time 3 :		Tomo time 7 :			
				Bucky RGDV								Tomo time 4 :		Tomo time 8 :			
				Tomo RGDV - switch related								Time setting for input at WA X21:18					
				Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready													
				switch related X11:3 --> Bucky - Tomo remote switchover RGDVs													

1) = has to be adjustet on site

[] = TDC

RGDV programming: example 2

BuckyDIAGNOST with Bucky-Controller
TH / TH2 / DigitalDiagnost

Name:

BuckyDIAGNOST with Bucky-Controller TH / TH2 / DigitalDiagnost														Name :							
Data Set A :														Bucky	Tomo	Bucky wall stand	Free cassette	RGDV5	RGDV6	RGDV7	RGDV8
Desk :														RGDV1	RGDV2	RGDV3	RGDV4				
- Room :														1	1	1	1				
- Tube :														1	1	1	1				
- Release circuit number :														do not care	do not care	do not care	do not care				
- Enable handswitch at generator desk :														yes	yes	yes	yes				
- Syncmaster present :														yes	yes	yes	yes				
- Exposure switch type :														double step	double step	double step	double step				
- Bucky format density correction :														0	0	0	0				
- Cone density correction :														0	0	0	0				
- Dose measurement input :														EZ X21	none / [EZ X21]	EZ X31	none				
- Dose measurement sensor :														Bucky amplimat	Bucky amplimat	Bucky amplimat	(Bucky amplimat)				
- Exposure series / Tomo movement :														no	no	no	no				
- Release delay :														enable	enable	enable	enable				
- Mounted radiographical controller :														Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost				
- Release circuit adaptation unit :														none	none	none	none				
- Mounted tomo extension :														none	none	none	none				
- Medium II format kV correction (dose equiv. steps) :														0	0	0	0				
- Medium II format density correction (6% steps) :														0	0	0	0				
- Medium II format mAs correction (6% steps) :														0	0	0	0				
- Small II format kV correction (dose equiv. steps) :														0	0	0	0				
- Small II format density correction (6% steps) :														0	0	0	0				
- Small II format mAs correction (6% steps) :														0	0	0	0				
Data Set B :																					
- Used for tomo :														no	yes	no	no				
- Used for fluoroscopy :														no	no	no	no				
- CT add on :														no	no	no	no				
- Disable time override :														no	no	no	no				
- Tube power factor :														100 %	100 %	100 %	100 %				
- kV steps :														Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)				
- mAs steps :														25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- mA steps :														25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- time steps :														25 % 1)	25 % 1)	25 % 1)	25 % 1)				
- Density steps :														12 % 1)	12 % 1)	12 % 1)	12 %				
- Density correction (6% steps) :														0	0	0	0				
- Underexposure display :														yes	yes	yes	yes				
- Tube overload protection :														on	on	on	on				

Bucky / Scopo 1WB / Decade Bucky 1 (WBX11):		Bucky / Tomo 1WA : Decade Bucky 1/2	WAX11	WAX12	Bucky / Tomo time
Bucky RGDV :	RGDV1 [] RGDV2 [] RGDV3 [] RGDV4 [] RGDV5 [] RGDV6 [] RGDV7 [] RGDV8 []	Tomo mode switch	----	----	Tomo time 1 : ----- Tomo time 2 : ----- Tomo time 3 : ----- Tomo time 4 : ----- Tomo time 5 : ----- Tomo time 6 : ----- Tomo time 7 : ----- Tomo time 8 : -----
		Tomo RGDV - switch related	----	----	Time setting for input at WA X21:18
		Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 ---> Bucky - Tomo remote switchover RGDV/s			

1) = has to be adjusted on site

$$[1] = \text{TDC}$$

Bucky/DIAGNOST TH any version with Bucky-Controller Generator equipped with / without decade adapt. unit WA Aux. for MCS (only) = RGDV4 combined with free cassette		Name :		Bucky		Tomo		Bucky wall stand		Free cassette		MCS			
Data Set A :		Desk :		RGDV1		RGDV2		RGDV3		RGDV4		RGDV4		RGDV6	
- Room :				1		1		1		1		1			
- Tube :				1		1		1		1		1			
- Release circuit number :				1		1		1		1		1			
- Enable handswitch at generator desk :				yes		yes		yes		yes		yes			
- Syncmaster present :				yes		yes		yes		yes		yes			
- Exposure switch type :				double step		double step		double step		double step		double step			
- Bucky format density correction :				0		0		0		0		0			
- Cone density correction :				0		0		0		0		0			
- Dose measurement input :				EZ X21		none / [EZ X21]		EZ X31		none		EZ X22			
- Dose measurement sensor :				Bucky amplimat		Bucky amplimat		Bucky amplimat		(Bucky amplimat)		Bucky amplimat			
- Exposure series / Tomo movement :				no		no		no		no		no			
- Release delay :				enable		enable		enable		enable		enable			
- Mounted radiographical controller :				Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost			
- Release circuit adaptation unit :				none		none		none		none		none			
- Mounted tomo extension :				none		none		none		none		none			
- Medium II format kV correction (dose equiv. steps) :				0		0		0		0		0			
- Medium II format density correction (6% steps) :				0		0		0		0		0			
- Medium II format mAs correction (6% steps) :				0		0		0		0		0			
- Small II format kV correction (dose equiv. steps) :				0		0		0		0		0			
- Small II format density correction (6% steps) :				0		0		0		0		0			
- Small II format mAs correction (6% steps) :				0		0		0		0		0			
Data Set B :															
- Used for tomo :				no		yes		no		no		no			
- Used for fluoroscopy :				no		no		no		no		no			
- CT add on :				no		no		no		no		no			
- Disable time override :				no		no		no		no		no			
- Tube power factor :				100 %		100 %		100 %		100 %		100 %			
- kV steps :				Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)		Dose equiv. 1)			
- mAs steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- mA steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- time steps :				25 % 1)		25 % 1)		25 % 1)		25 % 1)		25 % 1)			
- Density steps :				12 % 1)		12 % 1)		12 % 1)		12 %		12 % 1)			
- Density correction (6% steps) :				0		0		0		0		0			
- Underexposure display :				yes		yes		yes		yes		yes			
- Tube overload protection :				on		on		on		on		on			
Bucky / Scopo-1WB / Decade-Bucky 1 (WBX11) :		Bucky / Tomo-1WA : Decade-Bucky 1/2		WAX11		WAX12		Bucky / Tomo-time							
		Tomo mode switch		---		---		Tomo time 1 :		---		Tomo time 5 :		---	
		Bucky RGDV - switch related		---		---		Tomo time 2 :		---		Tomo time 6 :		---	
		Bucky RGDV		---		---		Tomo time 3 :		---		Tomo time 7 :		---	
		Bucky RGDV		---		---		Tomo time 4 :		---		Tomo time 8 :		---	
		Tomo RGDV - switch related		---		---		Time setting for input at WA X21:1							
		Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready													
		switch related X11:3 --> Bucky - Tomo remote switchover RGDVs													

1) = has to be adjustet on site

[] = TDC

RGDV programming: example 8

Name : Bucky/DIAGNOST TH any version Generator equipped with / without decade adapt. unit WA Aux. for MCS (only) = any of RGDV5... 8									
Data Set A :									
- Room :	RGDV1	Bucky	Tomo	Bucky wall stand	Free cassette	RGDV5	RGDV6	RGDV7	MCS
- Tube :	1		1	1	1				1
- Release circuit number :	1		1	1	1				1
- Enable handswitch at generator desk :	yes		yes	yes	yes				yes
- Syncmaster present :	yes		yes	yes	yes				no
- Exposure switch type :	double step		double step	double step	double step				double step
- Bucky format density correction :	0		0	0	0				0
- Cone density correction :	0		0	0	0				0
- Dose measurement input :	EZ X21		none / [EZ X21]	EZ X31	none				EZ X22
- Dose measurement sensor :	Bucky amplimat		Bucky amplimat	Bucky amplimat	(Bucky amplimat)				Bucky amplimat
- Exposure series / Tomo movement :	no		no	no	no				no
- Release delay :	enable		enable	enable	enable				enable
- Mounted radiographical controller :	Bucky contr. 1 / DigitalDiagnost		Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost	Bucky contr. 1 / DigitalDiagnost				none
- Release circuit adaptation unit :	none		none	none	none				none
- Mounted tomo extension :	none		none	none	none				none
- Medium II format kV correction (dose equiv. steps) :	0		0	0	0				0
- Medium II format density correction (6% steps) :	0		0	0	0				0
- Medium II format mAs correction (6% steps) :	0		0	0	0				0
- Small II format kV correction (dose equiv. steps) :	0		0	0	0				0
- Small II format density correction (6% steps) :	0		0	0	0				0
- Small II format mAs correction (6% steps) :	0		0	0	0				0
Data Set B :									
- Used for tomo :	no		yes	no	no				no
- Used for fluoroscopy :	no		no	no	no				no
- CT add on :	no		no	no	no				no
- Disable time override :	no		no	no	no				no
- Tube power factor :	100 %		100 %	100 %	100 %				100 %
- kV steps :	Dose equiv. 1)		Dose equiv. 1)	Dose equiv. 1)	Dose equiv. 1)				Dose equiv. 1)
- mAs steps :	25 % 1)		25 % 1)	25 % 1)	25 % 1)				25 % 1)
- mA steps :	25 % 1)		25 % 1)	25 % 1)	25 % 1)				25 % 1)
- time steps :	25 % 1)		25 % 1)	25 % 1)	25 % 1)				25 % 1)
- Density steps :	12 % 1)		12 % 1)	12 % 1)	12 % 1)				12 % 1)
- Density correction (6% steps) :	0		0	0	0				0
- Underexposure display :	yes		yes	yes	yes				yes
- Tube overload protection :	on		on	on	on				on
Bucky / Scopo 1WB / Decade Bucky 1 (WBX11) :-									
Bucky RGDV :	RGDV1 [] RGDV2 [] RGDV3 [] RGDV4 [] RGDV5 [] RGDV6 [] RGDV7 [] RGDV8 []		WAX11		WAX12		Bucky / Tomo time		
	Tomo mode switch		----		----		Tomo time 1 : -----		
	Bucky RGDV - switch related		----		----		Tomo time 2 : -----		
	Bucky RGDV		----		----		Tomo time 3 : -----		
	Bucky RGDV		----		----		Tomo time 4 : -----		
Tomo RGDV - switch related			----		----		Tomo time 5 : -----		
Tomo mode switch : X11:3 SL_XG_TO / Bucky RGDV : X11:1 Format + :10 Bucky ready / Tomo RGDV : X11:1 Format + :5 Tomo ready switch related X11:3 --> Bucky - Tomo remote switchover RGDVs			----		----		Tomo time 6 : -----		
			----		----		Tomo time 7 : -----		
			----		----		Tomo time 8 : -----		
			----		----		Tomo time 9 : -----		
			----		----		Tomo time 10 : -----		
			----		----		Tomo time 11 : -----		
			----		----		Tomo time 12 : -----		
			----		----		Tomo time 13 : -----		
			----		----		Tomo time 14 : -----		
			----		----		Tomo time 15 : -----		
			----		----		Tomo time 16 : -----		
			----		----		Tomo time 17 : -----		
			----		----		Tomo time 18 : -----		
			----		----		Tomo time 19 : -----		
			----		----		Tomo time 20 : -----		
			----		----		Tomo time 21 : -----		
			----		----		Tomo time 22 : -----		
			----		----		Tomo time 23 : -----		
			----		----		Tomo time 24 : -----		
			----		----		Tomo time 25 : -----		
			----		----		Tomo time 26 : -----		
			----		----		Tomo time 27 : -----		
			----		----		Tomo time 28 : -----		
			----		----		Tomo time 29 : -----		
			----		----		Tomo time 30 : -----		
			----		----		Tomo time 31 : -----		
			----		----		Tomo time 32 : -----		
			----		----		Tomo time 33 : -----		
			----		----		Tomo time 34 : -----		
			----		----		Tomo time 35 : -----		
			----		----		Tomo time 36 : -----		
			----		----		Tomo time 37 : -----		
			----		----		Tomo time 38 : -----		
			----		----		Tomo time 39 : -----		
			----		----		Tomo time 40 : -----		
			----		----		Tomo time 41 : -----		
			----		----		Tomo time 42 : -----		
			----		----		Tomo time 43 : -----		
			----		----		Tomo time 44 : -----		
			----		----		Tomo time 45 : -----		
			----		----		Tomo time 46 : -----		
			----		----		Tomo time 47 : -----		
			----		----		Tomo time 48 : -----		
			----		----		Tomo time 49 : -----		
			----		----		Tomo time 50 : -----		
			----		----		Tomo time 51 : -----		
			----		----		Tomo time 52 : -----		
			----		----		Tomo time 53 : -----		
			----		----		Tomo time 54 : -----		
			----		----		Tomo time 55 : -----		
			----		----		Tomo time 56 : -----		
			----		----		Tomo time 57 : -----		
			----		----		Tomo time 58 : -----		
			----		----		Tomo time 59 : -----		
			----		----		Tomo time 60 : -----		
			----		----		Tomo time 61 : -----		
			----		----		Tomo time 62 : -----		
			----		----		Tomo time 63 : -----		
			----		----		Tomo time 64 : -----		
			----		----		Tomo time 65 : -----		
			----		----		Tomo time 66 : -----		
			----		----		Tomo time 67 : -----		
			----		----		Tomo time 68 : -----		
			----		----		Tomo time 69 : -----		
			----		----		Tomo time 70 : -----		
			----		----		Tomo time 71 : -----		
			----		----		Tomo time 72 : -----		
			----		----		Tomo time 73 : -----		
			----		----		Tomo time 74 : -----		
			----		----		Tomo time 75 : -----		
			----		----		Tomo time 76 : -----		
			----		----		Tomo time 77 : -----		
			----		----		Tomo time 78 : -----		
			----		----		Tomo time 79 : -----		
			----		----		Tomo time 80 : -----		
			----		----		Tomo time 81 : -----		
			----		----		Tomo time 82 : -----		
			----		----		Tomo time 83 : -----		
			----		----		Tomo time 84 : -----		
			----		----		Tomo time 85 : -----		
			----		----		Tomo time 86 : -----		
			----		----		Tomo time 87 : -----		
			----		----		Tomo time 88 : -----		
			----		----		Tomo time 89 : -----		
			----		----		Tomo time 90 : -----		
			----		----		Tomo time 91 : -----		
			----		----		Tomo time 92 : -----		
			----		----		Tomo time 93 : -----		
			----		----		Tomo time 94 : -----		
			----		----		Tomo time 95 : -----		
			----		----		Tomo time 96 : -----		
			----		----		Tomo time 97 : -----		
			----		----		Tomo time 98 : -----		
			----		----		Tomo time 99 : -----		
			----		----		Tomo time 100 : -----		
			----		----		Tomo time 101 : -----		
			----		----		Tomo time 102 : -----		
			----		----		Tomo time 103 : -----		
			----		----		Tomo time 104 : -----		
			----		----		Tomo time 105 : -----		
			----		----		Tomo time 106 : -----		
			----		----		Tomo time 107 : -----		
			----		----		Tomo time 108 : -----		
			----		----		Tomo time 109 : -----		
			----		----		Tomo time 110 : -----		
			----		----		Tomo time 111 : -----		
			----		----		Tomo time 112 : -----		
			----		----		Tomo time 113 : -----		
			----		----		Tomo time 114 : -----		
			----		----		Tomo time 115 : -----		
			----		----		Tomo time 116 : -----		
			----		----		Tomo time 117 : -----		
			----		----		Tomo time 118 : -----		
			----		----		Tomo time 119 : -----		
			----		----		Tomo time 120 : -----		
			----		----		Tomo time 121 : -----		
			----		----		Tomo time 122 : -----		
			----		----		Tomo time 123 : -----		
			----		----		Tomo time 124 : -----		
			----		----		Tomo time 125 : -----		
			----		----		Tomo time 126 : -----		
			----		----		Tomo time 127 : -----		
			----		----		Tomo time 128 : -----		
			----		----		Tomo time 129 : -----		
			----		----		Tomo time 130 : -----		
			----		----		Tomo time 131 : -----		
			----		----		Tomo time 132 : -----		
			----		----		Tomo time 133 : -----		
			----		----		Tomo time 134 : -----		
			----		----		Tomo time 135 : -----		
			----		----		Tomo time 136 : -----		
			----		----		Tomo time 137 : -----		
			----		----		Tomo time 138 : -----		
			----		----		Tomo time 139 : -----		
			----		----		Tomo time 140 : -----		
			----		----		Tomo time 141 : -----		
			----		----		Tomo time 142 : -----		
			----		----		Tomo time 143 : -----		
			----		----		Tomo time 144 : -----		
			----		----		Tomo time 145 : -----		
			----		----		Tomo time 146 : -----		
			----		----		Tomo time 147 : -----		
			----		----		Tomo time 148 : -----		
			----		----		Tomo time 149 : -----		
			----		----		Tomo time 150 : -----		
			----		----		Tomo time 151 : -----		
			----		----		Tomo time 152 : -----		
			----		----		Tomo time 153 : -----		
			----		----		Tomo time 154 : -----		
			----		----		Tomo time 155 : -----		
			----		----		Tomo time 156 : -----		
			----		----		Tomo time 157 : -----		
			----		----		Tomo time 158 : -----		
			----		----		Tomo time 159 : -----		
			----		----		Tomo time 160 : -----		
			----		----		Tomo time 161 : -----		
			----		----		Tomo time 162 : -----		
			----		----		Tomo time 163 : -----		
			----		----		Tomo time 164 : -----		
			----		----		Tomo time 165 : -----		
			----		----		Tomo time 166 : -----		
			----		----		Tomo time 167 : -----		
			----		----		Tomo time 168 : -----		
			----		----		Tomo time 169 : -----		
			----		----		Tomo time 170 : -----		
			----		----		Tomo time 171 : -----		
			----		----		Tomo time 172 : -----		
			----		----		Tomo time 173 : -----		
			----		----		Tomo time 174 : -----		
			----		----		Tomo time 175 : -----		
			----		----		Tomo time 176 : -----		
			----		----		Tomo time 177 : -----		
			----		----		Tomo time 178 : -----		
			----		----		Tomo time 179 : -----		
			----		----		Tomo time 180 : -----		
			----		----		Tomo time 181 : -----		
			----		----		Tomo time 182 : -----		
			----		----		Tomo time 183 : -----		
			----		----		Tomo time 184 : -----		
			----		----		Tomo time 185 : -----		
			----		----		Tomo time 186 : -----		
			----		----		Tomo time 187 : -----		
			----		----		Tomo time 188 : -----		
			----		----		Tomo time 189 : -----		
			----		----		Tomo time 190 : -----		
			----		----		Tomo time 191 : -----		
			----		----		Tomo time 192 : -----		
			----		----		Tomo time 193 : -----		
			----		----		Tomo time 194 : -----		
			----		----		Tomo time 195 : -----		
			----		----		Tomo time 196 : -----		
			----		----		Tomo time 197 : -----		
			----		----		Tomo time 198 : -----		
			----		----		Tomo time 199 : -----		
			----		----		Tomo time 200 : -----		
			----		----		Tomo time 201 : -----		
			----		----		Tomo time 202 : -----		
			----		----		Tomo time 203 : -----		
			----		----		Tomo time 204 : -----		
			----		----		Tomo time 205 : -----		
			----		----		Tomo time 206 : -----		
			----		----		Tomo time 207 : -----		
			----		----		Tomo time 208 : -----		
			----		----		Tomo time 209 : -----		
			----		----		Tomo time 210 : -----		

1) = has to be adjustet on site

[] = TDC

RGDV programming: example 9